





• The most essential criteria for the long-term ECG recording is the quality of the electrode system. The custo tera program displays the current ECG of the patient on the PC screen directly after the recorder has been started. This guarantees a good ECG quality of the recording.

Recorder - custo flash 200

- Non-wearing microelectronics
- **OPTIONAL:** Recording of up to 3 independent ECG channels
- Quick data transmission with high-speed infrared technology
- Zero-line automatic system
- Automatic amplitude control
- Autostart
- Accumulation operation possible
- Size: 9.8 x 6.5 x 2.2 cm
- Weight incl. accumulator: 170 g

